

Product datasheet for RC223921

PTPN20A (NM_001042390) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PTPN20A (NM_001042390) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PTPN20A
Synonyms:	bA142I17.1; CT126; hPTPN20; PTPN20B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC223921 representing NM_001042390 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATTGTAACGATTATGAGGGAAATGACTCTGAAGCAGAAGACTTGAATTTTCAGGGAGACTTTGCCTT
CATCAAGTCAGGAAAACACACCTAGATCAAAGTTTTTGGAAAATAAAGTTAATTCAGAGAAGGTAAGTAACT
TTCTCTTCGGAATTTCCACATAATGATTATGAGGATGTTTTTGAAGAGCCTTCAGAAAGTGGCAGTGAT
CCCAGCATGTGGACAGCCAGAGGCCCTTCAGAAGAGACAGGTGGAGCAGTGAGGATGAGGAGGCTGCAG
GGCCATCACAGGCTCTCTCCCTCTACTTCTGATACGCGCAAATTTGTTTCTGAAGGAGAACTAGATCA
GTTGGCTCAGATTCGGCCATTAATATTCATTTTCATGAGCAGACGCCATCAAGGATTGTTTGGAAATC
CTTGAGGAAAAACAGCAGCGTATGATATCATGCAGGAATTTATGGCTTTAGAACTTAAGAATCTGCCTG
GTGAGTTCAACTCTGGGAATCAACCAAGCAACAGAGAAAAAACAGATACCGAGATATTCTTCCATATGA
TTCAACACGCGTTCTCTTGGAAAAAGCAAGGACTACATCAATGCTAGTTATATTAGAATAGTCAATTGT
GGAGAAGAGTATTTTTATCGCTACTCAAGGACCACTGCTGAGCACCATAGATGACTTTTGGCAAATGG
TGTTGGAAAATAATTCAAATGTTATTGCCATGATAACCAGAGAGATAGAAGGTGGAATTATCAAATGCTA
CCATACTGGCCATTTCTCTGAAGAAGCCATTGGAATTGAAACACTTCCGTGTATTCTGGAGAAGTAC
CAGATACTTCAATTTTCATCATTGCAATGTTTCAAGTTGTGGAGAAGTCCACGGAACTAGTCACTCTG
TAAAACAGTTGCAGTTCACCAAGTGGCCAGACCATGGCACTCCTGCCTCAGCAGATAGCTTCATAAAATA
TATTCGTTATGCAAGGAAGGCCACCTTACAGGACCCATGGTTGTTCACTGCAGTGCCGGCATAGGCCGG
ACAGGGGTGTTCTATGTGTGGATGTCGTGTTCTGTGCCATCGTAAAGAAGTGTTCATTCAACATCATGG
ATATAGTGGCCAAATGAGAGAACAACGTTCTGGCATGGTTCAAACGAAGGAGCAGTATCACTTTTGTTA
CGATATTGTGCTTGAAGTTCTTCGGAAACTTCTGACTTTGGAT

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC223921 representing NM_001042390
 Red=Cloning site Green=Tags(s)


MIVNDYEGNDSEAEDLNFRETLPSSSQENTPRSKVFENKVNSEKVKLSLRNFPHNDYEDVFEESPESGSD
 PSMWTARGPFRDRWSSSEDEEAAGPSQALSPLLSDTRKIVSEGELDQLAQIRPLIFNFHEQTAIKDCLKI
 LEEKTAAYDIMQEFMALELKNLPGFNSGNQPSNREKNRYRDILPYDSTRVPLGKSKDYINASYIRIVNC
 GEEYFYIATQGPLLSTIDDFWQMVLENNSNVIAMITREIEGGI IKCYHYWPI SLKKPLELKHFRVFL ENY
 QILQYFIIRMFQVVEKSTGTSHSVKQLQFTKWPDHGTPASADSF IKYIRYARKSHLTGPMVVHCSAGIGR
 TGVFLCVDVVFCAIVKNCSFNIMDIVAQMREQRSGMVQTKQYHFCYDIVLEVLRKLLTLD

SGP TRTRRLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8058_a11.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001042390

ORF Size: 1233 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001042390.1](#), [NP_001035849.1](#)

RefSeq Size: 2870 bp

RefSeq ORF: 1235 bp

Locus ID: 653129


Cytogenetics: 10q11.22

Protein Families: Druggable Genome

MW: 47.2 kDa

Gene Summary: The product of this gene belongs to the family of classical tyrosine-specific protein tyrosine phosphatases. Many protein tyrosine phosphatases have been shown to regulate fundamental cellular processes and several are mutated in human diseases. Chromosome 10q contains a segmental duplication resulting in multiple copies of the protein tyrosine phosphatase, non-receptor type 20 gene. The two nearly identical copies are designated as PTPN20A and PTPN20B. A third copy is only partially duplicated and contains a pseudogene, designated as PTPN20C. This gene encodes the more centromeric copy, PTPN20A. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC223921