

Product datasheet for **RC223923L3V**

Synaptotagmin 14 (SYT14) (NM_153262) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	Synaptotagmin 14
Synonyms:	SCAR11; sytXIV
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153262
ORF Size:	1665 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC223923).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_153262.1
RefSeq Size:	2960 bp
RefSeq ORF:	1668 bp
Locus ID:	255928
UniProt ID:	Q8NB59
Cytogenetics:	1q32.2



Domains: C2

Protein Families: Transmembrane

MW: 62.1 kDa

Gene Summary: This gene is a member of the synaptotagmin gene family and encodes a protein similar to other family members that mediate membrane trafficking in synaptic transmission. The encoded protein is a calcium-independent synaptotagmin. Mutations in this gene are a cause of autosomal recessive spinocerebellar ataxia-11 (SCAR11), and a t(1;3) translocation of this gene has been associated with neurodevelopmental abnormalities. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 4. [provided by RefSeq, Dec 2011]