

Product datasheet for **RC223978A1V**

Human ZNF713 (NM_182633) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ZNF713 (NM_182633) AAV Particle
Tag:	Myc-DDK
Symbol:	ZNF713
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC223978 representing NM_182633 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAGGAGGAAGAAATGAATGATGGCTCACAGATGGTGAGATCTCAGGAATCACTGACGTTTCAGGATG
TGGCCGTGGACTTCACCAGAGAGGAGTGGGACCAGCTGTACCCTGCCAAAAGAACCTCTATCGAGACGT
GATGCTGGAGAACTACAGGAATCTAGTTGCACTGGGTATCAGCTTTGTAAGCCAGAGGTAATCGCGCAG
TTGGAGCTAGAGGAAGAAATGGGTGATAGAAAGAGACAGCCTGCTGGATACTCATCCAGATGGAGAAAACA
GACCCGAAATCAAAAAGTCAACCACAAGCCAGAATATTCTGATGAAAATCAAACCCATGAGATGATAAT
GGAGAGACTCGCAGGAGACAGCTTCTGGTACTCCATCCTAGGAGGACTCTGGGATTTTGATTACCATCCA
GAGTTTAAACCAAGAAAACCAAGAGATATTTAGGACAAGTAACCTTGACCCACAAAAGATCACACAGG
AGAGAAGCCTTGAGTGAATAAATTTGCAGAAAATGTAATCTGAACTCAAACCTTATGCAGCAGAGAAT
TCCTTCCATTAATAATACCCCTGAATTCTGACACACAGGGAACAGCATCAAACATAATTCAGACTTGATT
TACTATCAGGGAAATTACGTAAGAGAGACTCCCTATGAATATAGTGAGTGTGGAAAAATCTTCAATCAAC
ATATTCTTCTTACTGATCATATTCATACTGCAGAGAAACCCAGTGAGTGTGGGAAGGCCCTTCAGCCACAC
CTCATCTCTTAGCCAGCCTCAGATGTTGCTTACAGGAGAGAAGCCCTATAAGTGTGATGAATGTGAAAAA
AGATTCAGCCAGAGGATACATCTCATTCAACATCAGAGAATTCACACAGGAGAAAAGCCTTTTATATGCA
ATGGATGTGGGAAGCCTTCGTCAGCATTATCCTTTTACTCAACATCTGAGGATTCATACTGGAGAAAA
GCCCTATAAATGAATCAATGTGGTAAAGCTTTTAGCCGCATCAGATCCCTTACTGAACATCATAGACTT
CATACCGGAGAGAAAACCTTACGAATGTGGTTTCTGTGGCAAAGCCTTCAGTCAGAGGACACATCTGAATC
AACATGAAAGAACTCATACAGGAGAGAAAACCTATAAATGTAATGAATGTGGGAAGCCTTTAGCCAGAG
TGCACACCTTAATCAACACAGGAAAATCCATACTCGGGAGAAATTATGTGAATATAAATGTGAGCAAAT
GTTCCGCCACAGTCCTTCATTTAGCAGCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC223978 representing NM_182633 Red=Cloning site Green=Tags(s) MEEEEENDGSQMVRSEQSLTFQDVAVDFTREEWDQLYPAQKNLYRDVMLENYRNLVALGYQLCKPEVIAQ LELEEEWVIERDSLDDTHPDGENRPEIKKSTTSQNI SDENQTHEMIMERLAGDSFWYSILGGLWDFDYHP EFNQENHKRYLGQVTLTHKKITQERSLECNKFAENCNLSNLMQQRIPSIKIPLNSDTQGNSIKHNSDLI YYQGNVRETPYEYSECGKIFNQHILLTDHIHTAEKPSECGKAFSHTSSLSQPQMLLTGEKPYKCDECGK RFSQRIHLIQHRIHTGEKPFICNGCGKA FRQHSSFTQHLRIHTGEKPYKCNQCGKAFSRITSLTEHHRL HTGEKPYECGFCGKAFSQRTHLNQHERHTHTGEKPYKCNECGKAFSQSAHLNQHRKIHTREKLCEYKCEQT VRHSPSFSST TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_182633
ORF Size:	1290 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_182633.1, NP_872439.1</u>
RefSeq Size:	2065 bp
RefSeq ORF:	1332 bp
Locus ID:	349075
Cytogenetics:	7p11.2
MW:	50 kDa