

## **Product datasheet for RC223991L4**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### MAGEA4 (NM\_001011550) Human Tagged Lenti ORF Clone

#### **Product data:**

**Product Type:** Expression Plasmids

Product Name: MAGEA4 (NM 001011550) Human Tagged Lenti ORF Clone

Tag: mGFP

Symbol: MAGEA4

Synonyms: CT1.4; MAGE-41; MAGE-X2; MAGE4; MAGE4A; MAGE4B

Mammalian Cell

Selection:

Puromycin

**Vector:** pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF in

Sequence:

The ORF insert of this clone is exactly the same as(RC223991).

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF

**ACCN:** NM\_001011550

ORF Size: 951 bp





#### MAGEA4 (NM\_001011550) Human Tagged Lenti ORF Clone - RC223991L4

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001011550.1</u>, <u>NP 001011550.1</u>

RefSeq Size: 1712 bp

RefSeg ORF: 954 bp

**Locus ID:** 4103

UniProt ID: P43358

Cytogenetics: Xq28

**MW:** 34.9 kDa


**Gene Summary:** This gene is a member of the MAGEA gene family. The members of this family encode

proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. Several variants encoding the

same protein have been found for this gene. [provided by RefSeq, Aug 2020]



# **Product images:**



Circular map for RC223991L4