

Product datasheet for **RC223991L4V**

MAGEA4 (NM_001011550) Human Tagged ORF Clone Lentiviral Particle

Product data:


Product Type:	Lentiviral Particles
Product Name:	MAGEA4 (NM_001011550) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAGEA4
Synonyms:	CT1.4; MAGE-41; MAGE-X2; MAGE4; MAGE4A; MAGE4B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001011550
ORF Size:	951 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC223991).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001011550.1 , NP_001011550.1
RefSeq Size:	1712 bp
RefSeq ORF:	954 bp
Locus ID:	4103
UniProt ID:	P43358
Cytogenetics:	Xq28
MW:	34.9 kDa



[View online »](#)

Gene Summary:

This gene is a member of the MAGEA gene family. The members of this family encode proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. Several variants encoding the same protein have been found for this gene. [provided by RefSeq, Aug 2020]

Product images:

[RC223991L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC223991L4V particle to overexpress human MAGEA4-mGFP fusion protein.