

Product datasheet for **RC223993A1V**

Human HAX1 (NM_001018837) AAV Particle

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: HAX1
Synonyms: HCLSBP1; HSIBP1; SCN3
Mammalian Cell None
Selection:
Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC223993 representing NM_001018837
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGAGCCTCTTTGATCTCTCCGGGGCTTTTCGGCTTTCCTGGACCTCGGAGCTTCAGCCAGGAGGAG
GGATACGTTTCCACGATAACTTCGGCTTTGATGACCTAGTACGAGATTTCAATAGCATCTTCAGCGATAT
GGGGGCTGGACCTTGCTTCCCATCCTCCTGAACTTCCAGGTCCTGAGTCAGAGACACCTGGTGAGAGA
CTACGGGAGGGACAGACACTTCGGGACTCAATGCTTAAGTATCCAGATAGTCACCAGCCAGGATCTTTG
GGGGGCTTGGAGAGTGATGAAGAAGTGAATCCCCCAACCAGCACCAGACTGGGGCTCCAGAGGCC
ATTCATAGGTTTGATGATGTATGGCCTATGGACCCCATCCTAGAACCCAGAGAGGACAATGATCTTGAT
TCCCAGGTTTCCAGGAGGGTCTTGGCCCGGTTCTACAGCCCAAGCCAAATCCTATTTCAGAGCATCT
CTGTGACCAAGATCACTAAACCAGATGGGATAGTGGAGGAGCGCCGGACTGTGGTGGACAGTGAGGGCCG
GACAGAGACTACAGTAACCCGACACGAAGCAGATAGCAGTCCTAGGGGTGATCCAGAATACCAAGACCT
CCAGCCCTGGATGATGCCTTTCCATCCTGGACTATTCTGGGACGTTGGTCCGGTCCCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC223993 representing NM_001018837
Red=Cloning site Green=Tags(s)

MSLFDLFRGFFGPGPRSFSPGGGIRFHDNFGFDDLVRDFNSIFSDMGAWTLPSHPPELPGPESETPGER
LREGQTLRDSMLKYPDSHQPRIFGGVLESDARSESPQAPDWGSQRPFHRFDDVWMPDPHPRTREDNDLD
SQVSQEGLGVLQPQPKSYFKSISVTKITKPDGIVEERRTVVDSEGRTEITVTRHEADSSPRGDPESPRP
PALDDAFSILDFLGRWFRSR

TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV**



Species:	Human
Serotype:	AAV-2
ACCN:	NM_001018837
ORF Size:	693 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001018837.1
RefSeq Size:	1052 bp
RefSeq ORF:	696 bp
Locus ID:	10456
UniProt ID:	O00165
Cytogenetics:	1q21.3
MW:	25.9 kDa