

Product datasheet for **RC224004**

SPACA5 (NM_205856) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SPACA5 (NM_205856) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SPACA5
Synonyms:	LYC5; LYZL5; PNP6288; SLLP-X; SLLP2; UNQ6288
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC224004 representing NM_205856 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGGCTGGGCACTGTGGTAGTGACCTTGGCCACGCTGATGGTTGTCCTGTGGATGCCAAGATCT
ATGAACGCTGCGAGCTGGCGCAAGACTGGAGAGAGCAGGGCTGAACGGCTACAAGGGCTACGGCGTTGG
AGACTGGCTGTGCATGGCTCATTATGAGAGTGGCTTTGACACCGCCTTCGTGGACCACAATCCTGATGGC
AGCAGTGAATATGGCATTTCCTCAACTGAATTCTGCCTGGTGGTGTGACAATGGCATTACACCCCAAGA
ACCTCTGCCACATGGATTGTCATGACCTGCTCAATCGCCATATTCTGGATGACATCAGGTGTGCCAAGCA
GATTGTGTCTCACAGAATGGGCTTTCTGCCTGGACTTCTGGAGGCTACACTGTTCTGGCCATGATTTA
TCTGAATGGCTCAAGGGTGTGATATGCATGTGAAAATTGATCCAAAATTCATCCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC224004 representing NM_205856 Red=Cloning site Green=Tags(s)
-------------------	--

MKAWGTVVVTLATLMVVTVDAKIYERCELAARLERAGLNGYKGYVGDWLCMAHYESGFDTAFVDHNPDG
SSEYGIQFLNSAWCDNGITPTKNLCHMDCHDLLNRHILDDIRCAKQIVSSQNL SAWTSWRHLHCSGHDL
SEWLKGCMDMHVKIDPKIHP

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Chromatograms:	https://cdn.origene.com/chromatograms/mk8044_a04.zip
----------------	---



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_205856

ORF Size: 477 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_205856.3](#), [NP_995328.2](#)


RefSeq Size: 732 bp

RefSeq ORF: 480 bp


Locus ID: 389852

UniProt ID: [Q96QH8](#)
 Cytogenetics: Xp11.23
 MW: 17.7 kDa

Product images:



Circular map for RC224004



Western blot validation of overexpression lysate (Cat# [LY404205]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC224004 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).