

Product datasheet for **RC224091A1V**

Human CD164L2 (NM_207397) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CD164L2 (NM_207397) AAV Particle
Tag: Myc-DDK
Symbol: CD164L2
Synonyms: EAPG6122; UNQ6122
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC224091 representing NM_207397
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGGCTCCGGGACCCGCGCCTTGCAGGACTGCGCTCTGTGGCGGCTGTTGCTGCCTCCTCTATGTG
CCCAGCTGGCTGTGGCTGGTAAAGGAGCTCGAGGCTTTGGGAGGGGAGCCCTGATCCGCCTGAATATCTG
GCCGGCGGTCCAAGGGGCTGCAACAGCTGGAGGTCTGTGAGCACTGCGTGGAGGGAGACAGAGCGCGC
AATCTCTCCAGCTGCATGTGGGAGCAGTGCCGGCCAGAGGAGCCAGGACACTGTGTGGCCCAATCTGAGG
TGGTCAAGGAAGGTTGCTCCATCTACAACCGCTCAGAGGCATGTCCAGCTGCTCACCACCACCCACCTA
TGAACCGAAGACAGTCAACACAGGGAGCCCCCAGTCCCTGAGGCCACAGCCCTGGATTGTACGGGGGCC
AGCTTTATCGGAGGTGTCGTGCTGTTGAGCCTACAGGCGGTGGCTTTCTTTGTGCTGCACTTCTCA
AGGCCAAGGACAGCACCTACCAGACGCTG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC224091 representing NM_207397
Red=Cloning site Green=Tags(s)

MEAPGPRALRTALCGGCCCLLLCAQLAVAGKGARGFGRGALIRLNIWPAVQGACKQLEVCEHCVEGDRAR
NLSSCMWEQCRPEEPGHCVAQSEVVKEGCSIYNRSEACPAHHHPTYPEKTVTTGSPVPPEAHSPGFDGA
SFIGGVVLVLSLQAVAFFVLHFLKAKDSTYQTL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_207397
ORF Size:	519 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_207397.1
RefSeq Size:	1517 bp
RefSeq ORF:	522 bp
Locus ID:	388611
UniProt ID:	Q6UWJ8
Cytogenetics:	1p36.11
MW:	18.2 kDa