

Product datasheet for **RC224091L1V**

CD164L2 (NM_207397) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CD164L2 (NM_207397) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CD164L2
Synonyms:	EAPG6122; UNQ6122
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_207397
ORF Size:	519 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224091).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_207397.1
RefSeq Size:	1517 bp
RefSeq ORF:	522 bp
Locus ID:	388611
UniProt ID:	Q6UWJ8
Cytogenetics:	1p36.11
Protein Families:	Transmembrane
MW:	18.2 kDa



[View online »](#)