

Product datasheet for **RC224136L3V**

SCUBE2 (NM_020974) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SCUBE2 (NM_020974) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SCUBE2
Synonyms:	CEGB1; CEGF1; CEGP1; scube/You
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020974
ORF Size:	2997 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224136).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020974.1 , NP_066025.1
RefSeq Size:	3737 bp
RefSeq ORF:	2916 bp
Locus ID:	57758
UniProt ID:	Q9NQ36
Cytogenetics:	11p15.4
Domains:	CUB, EGF_CA, EGF, EGF
Protein Families:	Transmembrane


[View online »](#)

MW: 106.6 kDa

Gene Summary: Lipid-binding protein required for SHH long-range signaling by binding to the dually lipid-modified SHH (ShhNp) and by promoting ShhNp mobilization, solubilization and release from the cell membrane (PubMed:22902404, PubMed:22677548). Acts by enhancing the proteolytic processing (shedding) of the lipid-modified N- and C- terminal of ShhNp at the cell surface (PubMed:24522195). Synergizes with DISP1 to increase SHH secretion (PubMed:22902404). Probable cell surface coreceptor for VEGFR2 involved in VEGFR2-mediated angiogenesis (PubMed:27834687).[UniProtKB/Swiss-Prot Function]