

Product datasheet for **RC224150**

FAM124A (NM_145019) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	FAM124A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



ORF Nucleotide Sequence: >RC224150 representing NM_145019
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGACCCAAAGGCGGGCGGCGGCGGAGGAGGACGACTGCGTGGACTCGGGCGCCGAGACCGGAGGT
 CCGACTACAGCCACCTGTCTCCACGAGCAATGGATTCTCCACCTGGCGGATGCCATGTTGCTGTCTCT
 TTAAGTACACACGGTGACAAAGAAAGGGAAGATGGAGCCTCCAGGTTGAACATACCTGGCCCTCAGGT
 GAGCTTTCCGTTGAAGAGGCGCAGGACCCCTTCTGGTCAGCATCCACATAATCGCAGACCCAGGGGAGT
 CCCAGCCCCGAGGAGGCCATCGACAACGTCCTGGCGTGGATCCACCCGACCTCCCGCTGTTCCGGGT
 GTCCGAGAGGCGGGCGTCCCGCGGCGGCGGAAGCCCCCAAGGGCGCTCAGCCAGCGCTGGCTGTGGTG
 CTGTTCTGCAGGAGGAGTACGGCGAAGAGCAGATCCTGCAGCTGCACCGCACACTGCAGCAGCCGCCCT
 GGCGCCACCACACCGAGCAGGTGCACGGCCGGTTCCTGCCCTACCTGCCCTGCAGCCAGGACTTCTT
 CACGCTGGCCCTGGGACGCCGCTTGGGCCATCCGGCCGTCGACTACGGCAAGGAAATCGTGCCTTC
 ACCGTCTACTGTGCTACGACAATATGCTGACAGCCTCAGGTTCTACAGCTGATTCTCCGAGGAACC
 CCAGCCAGAAGAAAGCGGACTTCTGCATCTTCCCTATTTTTTCAACCTGGATGTGGACATCCAGTTCTC
 CCTGAAAGACTGCCCTGTGACCACTGCCCGTGGCCACCGACTCCTCCGTGCTGGAGTCCGAGTGAGG
 GACATAGGCGAGCTCGTGCCTCTCTGCCAACCTTGACGCCCATCAGCGAGGGGCGCTGGCAGACGG
 AGGACCATGATGGGAACAAGATCCTCTACAGGCACAAAGGGTGCATAAGAAAGTTTCTAAACCTGGCAG
 AGTACATCATGCCTCCGAGAAGAAACGTCATTCCACTCCTTGGCGAGCACTGCTGTACCAAGCCATACA
 CCTGGCAGCAGCCAGCAGTCCCGCTCAACAGTCTCACCCGGGGCCATCCGACAGGCTGCCTCCTG
 GGCACAGCAGGAATTTGCCGACGAGCAACAGCACCCCAACCTCCCTGGTCTTCCAGAGAAGCAA
 GTCTTGTGTTTGTGTCACGGGAGGCCCTCCCTGGCTCCTCAGCTGAACCACAGTGGTTTTCAAAC
 ACAGGTGCCCCAGGGCACAGGGCATCAGAGTGGAGGCATGGCCACCTCCTCTCCATCGATGACCTAGAGG
 GGGCCAGGAGACAGCGTGGACACAGGCTGCGGTGTCTCATCGGACCTGTCTGTGGTCTCTGCATA
 TTCTGCACCCAGTAGGTTCTGCAGCAGTGGAGACACCCCTCCCTCCGAAAGATGCAGCAGCCACTGG
 GCAGCTCACAAAGGATTCAGGGAGGGACCACTGCCACTGTGACGAGGGTGACCACAGAGGCTCCTGGG
 CTTCCCTCCCTTTCTTACCAAAAGGTCTTCCAGCTCCTCAGCGACAGCTCGTGTGCTCCCCAGCTCC
 CAGCACCTCCACCCTCAGAGCTCCTCCACAGCTCCCATGCGATACCCCAAGTCAAGCAGACTGAT
 GGAGACATGCCACCACCCAGGCTGGCTGGCCCCGGGATAACGACATGGAGGAATCTACATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >RC224150 representing NM_145019
 Red=Cloning site Green=Tags(s)

MDPKAGGGGEEDDCVDSGAETGGSDYSHLSSTSNFSTWRMPCLFLYCTHGDKEKGRWSLQVEHTWPSG
 ELSVEEAQDPFLVSIHIIADPGESQPLQEAIDNVLAWIHPDLPLFRVSERRASRRRRKPPKGAQPALAVV
 LFLQEEYGEEQILQLHRTLQQPPWRHHHTEQVHGRFLPYLPCSQDFFTLAPGTPLWAIRPVHYGKEIVRF
 TVYCRYDNYADSLRFYQLILRRNPSQKKADFCIFPISNLDVDIQFSLKRLPCDQCPVPTDSSVLEFRVR
 DIGELVPLLPNPCSPISEGRWQTEDHDGNKILLQAQRVHKFKPKGRVHHASEKKRHSTPLPSTAVPSHT
 PGSSQSQPLNSPHPGPIRTGLPPGHQEQEFAGRANSTPNPPWSFQRSKSLFCLPTGGPSLASSAEPQWFSN
 TGAPGHRASEWRHGHLISIDDEGAQETDVTGLRLSSDLSVVSAYSAPSRFCSTVETPLPSERCSSHW
 AAHKDSREGPLPTVSRVTTEASWASLPFFTKRSSSSSATAAAPPAPSTSTLTDSSPQLPCDTPKVKQTD
 GDMPPPPGSAGPGDNDMEEFYI

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_145019

ORF Size: 1746 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:


1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_145019.2](#), [NP_659456.2](#)

RefSeq Size:	2104 bp
RefSeq ORF:	1749 bp
Locus ID:	220108
UniProt ID:	Q86V42
Cytogenetics:	13q14.3
MW:	64.2 kDa

Product images:



Circular map for RC224150



Western blot validation of overexpression lysate (Cat# [LY408081]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC224150 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).