

## Product datasheet for **RC224164L3V**

### SPOP (NM\_001007230) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | SPOP (NM_001007230) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | SPOP   |
| Synonyms:                 | BTBD32; NEDMACE; NEDMIDF; NSDVS1; NSDVS2; TEF2   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_001007230   |
| ORF Size:                 | 1122 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC224164).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001007230.1</a> , <a href="#">NP_001007231.1</a>  |
| RefSeq Size:              | 3034 bp  |
| RefSeq ORF:               | 1125 bp  |
| Locus ID:                 | 8405   |
| UniProt ID:               | <a href="#">O43791</a>   |
| Cytogenetics:             | 17q21.33   |
| MW:                       | 42.1 kDa   |



[View online »](#)

**Gene Summary:**

This gene encodes a protein that may modulate the transcriptional repression activities of death-associated protein 6 (DAXX), which interacts with histone deacetylase, core histones, and other histone-associated proteins. In mouse, the encoded protein binds to the putative leucine zipper domain of macroH2A1.2, a variant H2A histone that is enriched on inactivated X chromosomes. The BTB/POZ domain of this protein has been shown in other proteins to mediate transcriptional repression and to interact with components of histone deacetylase co-repressor complexes. Alternative splicing of this gene results in multiple transcript variants encoding the same protein. [provided by RefSeq, Jul 2008]