

Product datasheet for **RC224169L4V**

ZIK1 (NM_001010879) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ZIK1
Synonyms:	ZNF762
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001010879
ORF Size:	1461 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224169).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001010879.2
RefSeq Size:	2676 bp
RefSeq ORF:	1464 bp
Locus ID:	284307
UniProt ID:	Q3SY52
Cytogenetics:	19q13.43



MW: 54.8 kDa

Gene Summary: May be a transcriptional repressor.[UniProtKB/Swiss-Prot Function]