

Product datasheet for **RC224171L4V**

ANO8 (NM_020959) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	ANO8
Synonyms:	KIAA1623; TMEM16H
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020959
ORF Size:	3696 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224171).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_020959.1
RefSeq Size:	4150 bp
RefSeq ORF:	3699 bp
Locus ID:	57719
UniProt ID:	Q9HCE9
Cytogenetics:	19p13.11



Protein Families: Transmembrane

MW: 135.9 kDa

Gene Summary: Does not exhibit calcium-activated chloride channel (CaCC) activity.[UniProtKB/Swiss-Prot Function]