

Product datasheet for **RC224207L3V**

LCE2A (NM_178428) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	LCE2A
Synonyms:	LEP9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178428
ORF Size:	318 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224207).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_178428.3
RefSeq Size:	592 bp
RefSeq ORF:	321 bp
Locus ID:	353139
UniProt ID:	Q5TA79
Cytogenetics:	1q21.3



MW: 10.8 kDa

Gene Summary: Precursors of the cornified envelope of the stratum corneum.[UniProtKB/Swiss-Prot Function]