

## Product datasheet for **RC224209L2V**

### **ENOX2 (NM\_006375) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	ENOX2 (NM_006375) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ENOX2
Synonyms:	APK1; COVA1; tNOX
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006375
ORF Size:	1743 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224209).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006375.2</a> , <a href="#">NP_006366.2</a>
RefSeq Size:	4095 bp
RefSeq ORF:	1746 bp
Locus ID:	10495
UniProt ID:	<a href="#">Q16206</a>
Cytogenetics:	Xq26.1
Protein Families:	Secreted Protein
MW:	66.4 kDa



[View online »](#)

**Gene Summary:**

This gene is a tumor-specific member of the ECTO-NOX family of genes that encode cell surface NADH oxidases. The encoded protein has two enzymatic activities: catalysis of hydroquinone or NADH oxidation, and protein disulfide interchange. The protein also displays prion-like properties. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]