

Product datasheet for **RC224222L1V**

ENDOGL1 (EXOG) (NM_005107) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ENDOGL1 (EXOG) (NM_005107) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ENDOGL1
Synonyms:	ENDOGL1; ENDOGL2; ENGL; ENGL-a; ENGL-b; ENGLA; ENGLB
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005107
ORF Size:	1104 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224222).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005107.1
RefSeq Size:	1460 bp
RefSeq ORF:	1107 bp
Locus ID:	9941
UniProt ID:	Q9Y2C4
Cytogenetics:	3p22.2
Domains:	Endonuclease
Protein Families:	Druggable Genome, Transmembrane


[View online »](#)

Protein Pathways: Apoptosis

MW: 41.08 kDa

Gene Summary: This gene encodes an endo/exonuclease with 5'-3' exonuclease activity. The encoded enzyme catalyzes the hydrolysis of ester linkages at the 5' end of a nucleic acid chain. This enzyme is localized to the mitochondria and may play a role in programmed cell death. Alternatively spliced transcript variants have been described. A pseudogene exists on chromosome 18. [provided by RefSeq, Feb 2009]