

## Product datasheet for **RC224228L1V**

### SOX5 (NM\_006940) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SOX5 (NM_006940) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SOX5
Synonyms:	L-SOX5; L-SOX5B; L-SOX5F; LAMSHF
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_006940
ORF Size:	2289 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224228).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006940.4</a>
RefSeq Size:	4333 bp
RefSeq ORF:	2292 bp
Locus ID:	6660
UniProt ID:	<a href="#">P35711</a>
Cytogenetics:	12p12.1
Domains:	HMG
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 83.8 kDa

**Gene Summary:** This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. The encoded protein may play a role in chondrogenesis. A pseudogene of this gene is located on chromosome 8. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]