

Product datasheet for **RC224229A1V**

Human WWP2 (NM_199423) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human WWP2 (NM_199423) AAV Particle
Tag:	Myc-DDK
Symbol:	WWP2
Synonyms:	AIP2; WWp2-like
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC224229 representing NM_199423 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCATCTGCCAGCTCTAGCCGGGCAGGAGTGGCCCTGCCTTTTGAGAAGTCTCAGCTCACTTTGAAAG
TGGTGTCCGCAAAGCCCAAGGTGCATAATCGTCAACCTCGAATTAACCTACGTGGAGGTGGCGGTGGA
TGGACTCCCCAGTGAGACCAAGAAGACTGGGAAGCGCATTGGGAGCTCTGAGCTTCTCTGAATGAGATC
ATCATTGTAATGTCACGGCACAGAGTCATTTAGATTTAAAGGTCTGGAGCTGCCATACCTTGAGAAATG
AACTGCTAGGCACCGCATCTGTCAACCTCTCCAACGCTTGAAGAACAATGGGGGCAAAATGGAGAACAT
GCAGCTGACCCTGAACCTGCAGACGGAGAACAAGGCAGCGTTGTCTCAGGCGGAGAGCTGACAATTTTC
CTGGACGGGCCAACTGTTGATCTGGGAAATGTGCCTAATGGCAGTGCCCTGACAGATGGATCACAGCTGC
CTTCGAGAGACTCCAGTGGAACAGCAGTAGCTCCAGAGAACCGGCACCGCCCCAGCACAACTGCTT
TGGTGGAAGATCCCGGACGCACAGACATTGGGTGCTTCAGCCAGAACAACCCAGCAACCGGCGAGCAA
AGCCCCGGTGCTCGGAGCCGGCACCGCCAGCCGTCAGAACTCAGGCCACAGTGGCTTGCCCAATGGCA
CAGTGAATGATGAACCCACAACAGCCACTGATCCCGAAGAACCTCCGTTGTTGGTGTGACGTCCCCACC
TGCTGCACCCTTGAGTGTGACCCGAATCCCAACAGACTTCTCTCCCTGCCCCAGCCACACCGGCTGAA
GGAGAGGAACCCAGCACTTCGGGTACACAGCAGCTCCAGCGGCTGCCAGGCCCGGACGCTCTGCCTG
CTGGATGGGAACAGCGAGAGCTGCCAACGGACGTGTCTATTATGTTGACCACAATACCAAGACCACCAC
CTGGGAGCGGCCCTTCTCCAGGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC224229 representing NM_199423 Red=Cloning site Green=Tags(s) MASASSSRAGVALPFEKSQLTLKVVSAPKPKVHNRQPRINSYVEVAVDGLPSETKKTKGRIGSSELLWNEI IILNVT AQSHLDLKVWSCHTLRNELLGTASVNL SNVLKNNGGKMENMQLTLNLQTENKGSVVSGGELTIF LDGPTVDLGNVPNGSALTDGSQLPSRDSSGTAVAPENRHQPPSTNCFGGRSRTHRHSGASARTTPATGEQ SPGARSRRHQPVKNSGHSGLANGTVNDEPTTATDPEEPSVVGVTSPPAAPLSVTPNPNTTSLPAPATPAE GEEPSTSGTQQLPAAAQAPDALPAGWEQRELPNGRVYYVDHNTKTTTWERPLPPG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_199423
ORF Size:	1005 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_199423.1, NP_955455.1</u>
RefSeq Size:	2659 bp
RefSeq ORF:	1007 bp
Locus ID:	11060
Cytogenetics:	16q22.1
MW:	35.1 kDa