

Product datasheet for **RC224258L3V**

PABPN1L (NM_001080487) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PABPN1L (NM_001080487) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PABPN1L
Synonyms:	ePABP2; PABPNL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001080487
ORF Size:	834 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224258).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001080487.1
RefSeq Size:	870 bp
RefSeq ORF:	837 bp
Locus ID:	390748
UniProt ID:	A6NDY0
Cytogenetics:	16q24.3
MW:	31.3 kDa
Gene Summary:	Binds the poly(A) tail of mRNA.[UniProtKB/Swiss-Prot Function]



[View online »](#)