

Product datasheet for **RC224259L4V**

NTN3 (NM_006181) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	NTN3
Synonyms:	NTN2L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006181
ORF Size:	1740 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224259).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_006181.1
RefSeq Size:	1986 bp
RefSeq ORF:	1743 bp
Locus ID:	4917
UniProt ID:	O00634
Cytogenetics:	16p13.3



Protein Pathways: Axon guidance

MW: 61.3 kDa

Gene Summary: Netrins control guidance of CNS commissural axons and peripheral motor axons.
[UniProtKB/Swiss-Prot Function]