

Product datasheet for **RC224260A1V**

Human BEND6 (NM_152731) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human BEND6 (NM_152731) AAV Particle
Tag:	Myc-DDK
Symbol:	BEND6
Synonyms:	C6orf65
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC224260 representing NM_152731 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCAGAAGATCGTGCAGACAGATGAAATTACCAATACACAAGCTTTAGAAAAGGAAAGAGGAAAAGAA
CAGAGACAATGGACTCAGAAAATGCAATAGTGACATGGATAAAGGACAGAGAGACCCATATTCGGGAAA
TGCCTTTCTGCCTGGTGAAGCTCCAGTGAGGATGAAGAGCCTTTAGCAGAATTGTCAAAGGAAGAATTG
TGCGCCAAAATAAAAGCCTGAAAGAAAACTAACAAACACCCGAAAGAAAACAGCCGACTTCGACAGT
CTTTGGTCATGCTTCAAGTGTTACCACAAGCAGTCACCCAGTTTGAAGAATTGGTTGGTATGGCCGAGGC
TCTGCTTAAGGGTGGGGGAACCATGTCTACATCTGCATCCACCCTCTGGAGAGCAACAAACAACCTCTCG
CCAGATTCATTTGCCTCAACATGCAGTAATTCTAATTCTAACTCCAGTTCACCCAGTTTCCTTAAAGCCTG
AGGAAGAGCATCAGACTGATGAGAAACAGTTCAGATTGAAAAATGGCAGATTGCCGTTGTAACAAGAG
CAAGCCTCAGAAGTTTATTAATGATTTAATGCAAGTACTTTACACAAATGAATACATGGCCACTCACAGC
CTGACAGGGGCAAAATCCTCTACTTCAAGGGACAAAGCTGTAAAACCAGCTATGAATCAGAATGAAGTCC
AAGAAATCATAGGAGTAACAAAACAATTATTTCCCAATACGGATGATGTTTCAATTAGGAGAATGATAGG
GCAAAAGCTAAACAACGTACCAAGAAGCCAAATTTAAGCAAAAATCTTAACCTCTCAGGATATTTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC224260 representing NM_152731 Red=Cloning site Green=Tags(s) MQKIVQTDEITNTQAFRKGKRKRTETMDSENANSMDKGQRPYSGNAFLPGESSSEDEEPLAELSKEEL CAKIKSLKEKL TNTRKENSRLRQSLVMLQVLPQAVTQFEELVGMAEALLKGGGTMTSASTLWRATNNSS PDSFASTCSNSNSNSSSPVSLKPEEEHQTDKQFQIEKWQIARCNKSKPQKFINDLMQVLYTNEYMATHS LTGAKSSTSRDKAVKPAMNQNEVQEIIIGVTKQLFPNTDDVSIRRMIGQKLNCTKKPNLSKNLNSQDIK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_152731
ORF Size:	837 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_152731.2</u>
RefSeq Size:	2727 bp
RefSeq ORF:	840 bp
Locus ID:	221336
UniProt ID:	<u>Q5SZJ8</u>
Cytogenetics:	6p12.1
MW:	31 kDa