

Product datasheet for **RC224266L4V**

TMTC3 (NM_181783) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	TMTC3
Synonyms:	LIS8; SMILE
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_181783
ORF Size:	2742 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224266).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_181783.1
RefSeq Size:	7203 bp
RefSeq ORF:	2745 bp
Locus ID:	160418
UniProt ID:	Q6ZXV5
Cytogenetics:	12q21.32



Protein Families: Transmembrane

MW: 103.9 kDa

Gene Summary: This gene encodes a protein that belongs to the transmembrane and tetratricopeptide repeat-containing protein family. [provided by RefSeq, May 2010]