

## Product datasheet for **RC224267A1V**

### Human TMEM18 (NM\_152834) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human TMEM18 (NM\_152834) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** TMEM18  
**Synonyms:** IncND  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC224267 representing NM\_152834  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCCGTCGCGCTTCTCTGTCTGCTCAGCTCTTTCCCGTCAGCATCCCAGCCGTGCTCACGCAGACGGACTGGA  
CTGAGCCCTGGCTCATGGGGCTGGCCACCTCCACGCGCTCTGCGTGCTCCTCACCTGCTTGCTCCCG  
AAGCTACAGACTACAGATCGGGCACTTTCTGTGTCTAGTCATCTTAGTCTACTGTGCTGAATACATCAAT  
GAGGCGGCTGCGATGAAGTGGAGATTATTTTCGAAATACCAGTATTTTCGACTCCAGGGGGATGTTTCATTT  
CTATAGTATTTTCAGCCCACTGCTGGTGAATGCCATGATCATTGTGGTTATGTGGGTATGGAAGACTTT  
GAATGTGATGACTGACCTGAAGAATGCACAAGAGAGAAGAAAGAAAAGAAAAGGAGAAGGAAAGAAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC224267 representing NM\_152834  
Red=Cloning site Green=Tags(s)

MPSAFSVSSFPVSI PAVL TQTDWTEPWL MGLATFHALCVLLTCLSSRSYRLQIGHFLCLVILVYCAEYIN  
EAAAMNWR LFSKYQYFDSRGMFISIVFSAPLLVNAMIIIVMWWKTLNVMTDLKNAQERRRKEKRRRRKED

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_152834



[View online »](#)

<b>ORF Size:</b>	420 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_152834.2</a></u>
<b>RefSeq Size:</b>	2138 bp
<b>RefSeq ORF:</b>	423 bp
<b>Locus ID:</b>	129787
<b>UniProt ID:</b>	<u><a href="#">Q96B42</a></u>
<b>Cytogenetics:</b>	2p25.3
<b>MW:</b>	16.1 kDa