

Product datasheet for **RC224286**

BAGE5 (NM_182484) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	BAGE5
Synonyms:	CT2.5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC224286 representing NM_182484
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAGCCGGAGCGGTTTTCTGGCATTGTCTGCCAGCTGCTCCAAGCCAGACTGATGAAGGAGGAGT
CCCCTGTGGTGAGCTGGAGTTGGAGCCTGAAGACGGCACAGCTCTGTGCTTCATCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC224286 representing NM_182484
Red=Cloning site Green=Tags(s)

MAAGAVFLALSAQLLQARLMKEESPVSWRLEPEDGTALCFIF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8043_f07.zip

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_182484

ORF Size: 129 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)


OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_182484.2</u>
RefSeq Size:	1840 bp
RefSeq ORF:	132 bp
Locus ID:	85316
UniProt ID:	<u>Q86Y27</u>
Cytogenetics:	13cen
MW:	5.2 kDa
Gene Summary:	Unknown. Candidate gene encoding tumor antigens.[UniProtKB/Swiss-Prot Function]

Product images:

Circular map for RC224286