

Product datasheet for **RC224301**

KRTAP4-11 (NM_033059) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	KRTAP4-11
Synonyms:	KAP4.11; KAP4.14; KRTAP4-14; KRTAP4.14
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC224301 representing NM_033059
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGTAAACTCCTGTTGTGGCTCCGTGTGCTCTACCAAGGCTGTGGCCAAGACCTCTGCCAGGAGACCT
GCTGCCGCCCCAGCTGCTGTGAGACCACCTGCTGCAGGACCACCTACTGTGCCCCAGCTGCTGTGTGTC
CAGCTGCTGCAGGCCCCAGTGTGCCAGTCTGTGTGCTGCCAGCCACCTGCTGCCGCCCCAGATGCTGC
ATCTCCAGCTGCTGTGCCCCAGTGTGTGTGCCAGCTGCTGCAAGCCCCAGTGTGCCAGTCTATGT
GCTGCCAGCCCACTTGTGCCGCCCCAGATGCTGCATCTCCAGCTGCTGTGCCCCAGCTGCTGTGTGTC
CAGCTGCTGCAGACCCAGTGTGCCAGTCTGTGTGCTGCCAGCCACCTGCTGCCACCCAGCTGCAGC
ATCTCCAGCTGCTGCCGCCCCCTTTGCTGTGAATCCAGCTGCTGCCGCCCCCTGCTGCTGCCTGCCAG
TCTGTGGCGGAGTCTCCTGCCACCACTTGTATCGCCCAACCTGTGTATCTCCAGCTGCCCCGCC
CTTGTGCTGTGCCTCCTTGTGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC224301 representing NM_033059
Red=Cloning site Green=Tags(s)

MVNSCCGSVCSHQGCGQDLQETCCRPSCCETTCCRTTYCRPSCCVSSCCRPQCCQSVCCQPTCCRPCC
ISSCCRPSCCVSSCKPQCCQSMCCQPTCCRPCCISSCCRPSCCVSSCCRPQCCQSVCCQPTCCHPSCS
ISSCCRPSCCESSCCRPCCCLRPVCGGVSCHTTCYRPTCVISSCPRLCCASSCC


TRTRPLEQKLI SEEDLAANDILDYKDDDDKV



Chromatograms: https://cdn.origene.com/chromatograms/mk8043_e03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_033059

ORF Size: 585 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_033059.2, NP_149048.1</u>
RefSeq Size:	1219 bp
RefSeq ORF:	588 bp
Locus ID:	653240
UniProt ID:	<u>Q9BYQ6</u>
Cytogenetics:	17q21.2
MW:	20.7 kDa
Gene Summary:	This protein is a member of the keratin-associated protein (KAP) family. The KAP proteins form a matrix of keratin intermediate filaments which contribute to the structure of hair fibers. KAP family members appear to have unique, family-specific amino- and carboxyl-terminal regions and are subdivided into three multi-gene families according to amino acid composition: the high sulfur, the ultrahigh sulfur, and the high tyrosine/glycine KAPs. This protein is a member of the ultrahigh sulfur KAP family and the gene is localized to a cluster of KAPs at 17q12-q21. [provided by RefSeq, Mar 2009]

Product images:



Circular map for RC224301