

Product datasheet for **RC224309L4V**

C14orf181 (NM_207442) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	C14orf181
Synonyms:	FLJ39779
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_207442
ORF Size:	498 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224309).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_207442.1 , NP_997325.2
RefSeq Size:	1727 bp
RefSeq ORF:	500 bp
Locus ID:	400223
Cytogenetics:	14q24.1
MW:	17.1 kDa



Gene Summary:

DISCONTINUED: This record has been withdrawn by NCBI and HGNC. This record was defined by AK097098.1, which is not sufficient evidence to define a distinct locus.