

Product datasheet for **RC224314A1V**

Human **SSX4** (NM_005636) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SSX4 (NM_005636) AAV Particle
Tag: Myc-DDK
Symbol: SSX4
Synonyms: CT5.4
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC224314 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAACGGAGACGACGCTTTGCAAGGAGACCCAGGGATGATGCTCAAATATCAGAGAAGTTACGAAAGG
CCTTCGATGATATTGCCAAATACTTCTCTAAGAAAGAGTGGGAAAAGATGAAATCCTCGGAGAAAATCGT
CTATGTGTATATGAAGCTAACTATGAGGTCATGACTAACTAGGTTTCAAGGTCACCTCCCACCTTTC
ATGCGTAGTAAACGGGCTGCAGACTTCCACGGGAATGATTTTGGTAACGATCGAAACACAGGAATCAGG
TTGAACGTCCTCAGATGACTTTCGGCAGCCTCCAGAGAATCTTCCGAAGATCATGCCAAGAAGCCAGC
AGAGGAAGAAAATGGTTTGAAGGAAGTGCCAGAGGCATCTGGCCACAAAATGATGGGAAACAGCTGTGC
CCCCGGGAAATCCAAGTACCTTGGAGAAGATCAACAAGACATCTGGACCCAAAAGGGGAAACATGCCT
GGACCCACAGACTGCGTGAGAGAAAGCAGCTGGTGGTTTATGAAGAGATCAGCGACCCTGAGGAAGATGA
CGAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC224314 protein sequence
Red=Cloning site Green=Tags(s)

MNGDDAFARRPRDDAQISEKLKAFDDIAKYFSKKEWEKMSSEKIVVYVMKLNVEVMTKLGFKVTLPF
MRSKRAADFHGNDFGNDRNHRNQVERPQMTFGSLQRIFFPKIMPKKPAEEENGLKEVPEASGPQNDGKQLC
PPGNPSTLEKINKTSGPKRGKHAWTHRLRERKQLVVYEEISDPEEDDE

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_005636
ORF Size:	564 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_005636.1
RefSeq Size:	1250 bp
RefSeq ORF:	567 bp
Locus ID:	6759
UniProt ID:	O60224
Cytogenetics:	Xp11.23
MW:	21.9 kDa