

Product datasheet for **RC224320L4V**

RIMS4 (NM_182970) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RIMS4 (NM_182970) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RIMS4
Synonyms:	C20orf190; RIM-4; RIM 4; RIM4; RIM4-gamma; RIM4gamma
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_182970
ORF Size:	807 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224320).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_182970.2
RefSeq Size:	5136 bp
RefSeq ORF:	810 bp
Locus ID:	140730
UniProt ID:	Q9H426
Cytogenetics:	20q13.12
MW:	29.1 kDa
Gene Summary:	Regulates synaptic membrane exocytosis.[UniProtKB/Swiss-Prot Function]


[View online »](#)