

Product datasheet for **RC224334L4V**

BAGE4 (NM_181704) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	BAGE4
Synonyms:	CT2.4; MLL3P
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_181704
ORF Size:	117 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224334).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_181704.1 , NP_859055.1
RefSeq Size:	1840 bp
RefSeq ORF:	120 bp
Locus ID:	85317
UniProt ID:	Q86Y28
Cytogenetics:	2p11.1



Protein Families: Druggable Genome

MW: 4.7 kDa

Gene Summary: Unknown. Candidate gene encoding tumor antigens.[UniProtKB/Swiss-Prot Function]