

## Product datasheet for **RC224360**

### **C2orf88 (NM\_001042521) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** C2orf88 (NM\_001042521) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** C2orf88  
**Synonyms:** smAKAP  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC224360 ORF sequence  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGCTGCATGAAATCAAAGCAAACCTTCCCATTTCTACCATATATGAAGGTGAGAAGCAGCATGAGA  
GTGAAGAACCCTTTATGCCAGAAGAGAGATGTCTACCTAGGATGGCTTCTCCAGTTAATGTCAAAGAGGA  
AGTGAAGGAACCTCCAGGGACCAATATTGTGATCTTGAATATGCACACCGCCTGTCTCAGGATATCTTG  
TGTGATGCCTTGACGAATGGGCATGCAATAACATCAAGTACCATGACATTCCATACATTGAGAGTGAGG  
GGCCT

**ACGGCT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC224360 protein sequence  
**Red**=Cloning site **Green**=Tags(s)  
MGMKSKQTFPFPTIYEGEKQHESEEPFMPEERCLPRMASPVNVKKEEVEKPPGTNIVILEYAHRLSQDIL  
CDALQQWACNNIKYHDIPYIESEGP

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk6047\\_a03.zip](https://cdn.origene.com/chromatograms/mk6047_a03.zip)

**Restriction Sites:** SgfI-MluI



[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001042521

**ORF Size:** 285 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001042521.1](#), [NP\\_001035986.1](#)

**RefSeq Size:** 4019 bp

**RefSeq ORF:** 288 bp

**Locus ID:** 84281


UniProt ID: [Q9BSF0](#)

Cytogenetics: 2q32.2


MW: 11 kDa

Gene Summary: Binds to type I regulatory subunits of protein kinase A (PKA-R1) and may anchor/target them to the plasma membrane.[UniProtKB/Swiss-Prot Function]


**Product images:**



Circular map for RC224360



Western blot validation of overexpression lysate (Cat# [LY420959]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC224360 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified C2orf88 protein (Cat# [TP324360]). The protein was produced from HEK293T cells transfected with C2orf88 cDNA clone (Cat# RC224360) using MegaTran 2.0 (Cat# [TT210002]).