

Product datasheet for **RC224364A1V**

Human PRRG2 (NM_000951) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PRRG2 (NM_000951) AAV Particle
Tag: Myc-DDK
Symbol: PRRG2
Synonyms: PRGP2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC224364 representing NM_000951
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGGGGCCACCCCTCTCTGCTGCTGCTATATATGGCATTAAACCACCTGCCTGGATACTTCACCCAGTG
AGGAGACAGACCAAGAAGTCTTCCTGGGTCCCCAGAGGCCAGAGCTTCCTGAGTAGCCATACCCGGAT
TCCAAGAGCCAACCACTGGGACCTGGAGCTGCTCACACCAGGGAACCTGGAACGGGAGTGCTGGAAGAG
AGGTGTTCTCTGGGAAGAGGCCAGGGAGTATTTTGAGGACAACACTCTCACGGAGCGCTTTTGGGAGAGCT
ACATCTACAATGGCAAAGGAGGGCGTGGACGAGTGGATGTGGCCAGCCTGGCTGTGGGGCTGACAGGTGG
CATCCTGCTCATTGTCTGGCCGGCCTGGGAGCCTTTTGGTATCTGCGCTGGCGACAGCACCAGAGGCCAG
CAGCCCTGTCCCAAGAGGCCGGGCTCATTAGCCCTCTGAGTCCTTTGAACCTCTGGGCCACCCAGCGC
CCCTGCCTCCACCCCAACCCCAACCCCAAGGCCTCCACCTATGAGCAGGCGCTGGCAGCCTCTGGGGT
ACACGACGCACCTCCACCCCTACACAGCCTCAGGAGGCCTCAC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC224364 representing NM_000951
Red=Cloning site Green=Tags(s)

MRGHPSLLLLLYMALTTCLDTSPSEETDQEVFLGPPEAQSFSSHTRIPRANHWDLELLTPGNLERECLEE
RCSWEEAREYFEDNTLTERFWESIYNGKGRGRVDVASLAVGLTGGILLIVLAGLAFWYLRWRQHRGQ
QPCPQEAGLISPLSPLNPLGPPTPLPPPPPPPLPTYEQALAASGVHDAPPPYTSLRPH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_000951
ORF Size:	606 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_000951.1</u>
RefSeq Size:	1167 bp
RefSeq ORF:	609 bp
Locus ID:	5639
UniProt ID:	<u>O14669</u>
Cytogenetics:	19q13.33
MW:	22.2 kDa