

Product datasheet for **RC224365L4V**

PAXBP1 (NM_016631) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PAXBP1 (NM_016631) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PAXBP1
Synonyms:	BM020; C21orf66; FSAP105; GCFC; GCFC1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016631
ORF Size:	2751 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224365).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016631.2
RefSeq Size:	2903 bp
RefSeq ORF:	2754 bp
Locus ID:	94104
UniProt ID:	Q9Y5B6
Cytogenetics:	21q22.11
Protein Families:	Transcription Factors
MW:	104.6 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein that may bind to GC-rich DNA sequences, which suggests its involvement in the regulation of transcription. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jun 2009]