

Product datasheet for **RC224386L4V**

CWH43 (NM_025087) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	CWH43
Synonyms:	CWH43-C; PGAP2IP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025087
ORF Size:	2097 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224386).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_025087.1
RefSeq Size:	2486 bp
RefSeq ORF:	2100 bp
Locus ID:	80157
UniProt ID:	Q9H720
Cytogenetics:	4p11



Protein Families:	Transmembrane
MW:	78.4 kDa
Gene Summary:	Involved in lipid remodeling during GPI-anchor maturation.[UniProtKB/Swiss-Prot Function]