

Product datasheet for RC224406L3V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

VEGF Receptor 1 (FLT1) (NM 002019) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: VEGF Receptor 1 (FLT1) (NM_002019) Human Tagged ORF Clone Lentiviral Particle

Symbol: FLT1

Synonyms: FLT; FLT-1; VEGFR-1; VEGFR1

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK ACCN: NM_002019

ORF Size: 4014 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC224406).

Sequence:

Cytogenetics:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through

naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 002019.2, NP 002010.2

 RefSeq Size:
 5777 bp

 RefSeq ORF:
 4017 bp

 Locus ID:
 2321

 UniProt ID:
 P17948

Domains: pkinase, TyrKc, S_TKc, ig, IGv, IGc2, IG

13q12.3

Protein Families: Druggable Genome, Protein Kinase, Secreted Protein





Protein Pathways: Cytokine-cytokine receptor interaction, Endocytosis, Focal adhesion

MW: 150.77 kDa

Gene Summary: This gene encodes a member of the vascular endothelial growth factor receptor (VEGFR)

family. VEGFR family members are receptor tyrosine kinases (RTKs) which contain an extracellular ligand-binding region with seven immunoglobulin (Ig)-like domains, a transmembrane segment, and a tyrosine kinase (TK) domain within the cytoplasmic domain. This protein binds to VEGFR-A, VEGFR-B and placental growth factor and plays an important role in angiogenesis and vasculogenesis. Expression of this receptor is found in vascular endothelial cells, placental trophoblast cells and peripheral blood monocytes. Multiple transcript variants encoding different isoforms have been found for this gene. Isoforms include a full-length transmembrane receptor isoform and shortened, soluble isoforms. The

soluble isoforms are associated with the onset of pre-eclampsia.[provided by RefSeq, May

2009]