

Product datasheet for **RC224407**

N42L1 (N4BP2L1) (NM_052818) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	N42L1
Synonyms:	CG018
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC224407 representing NM_052818
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAGGACAGTTTCCTTCAATCTTTGGGAGGCTGAGCCTCCAGCCCCAGCAGCAGCAGCAGCGGCAGC
GGCCGCCCCGGCCGCCCGCGGGGACACCTCCTCGCCGCCACAGCTTTAGGAAACACCTCTACCTCCT
GCGAGGCCTCCCGGGCTCCGGGAAACTACACTGGCCAGACAATTGCAGCATGACTTTCCAGGGCCCTG
ATTTTCAGCACGGATGATTTTTCTTCAGGGAAGATGGTGCCTATGAGTTCAATCCTGACTTCTCTGGAGG
AAGCTCATGAATGGAACCAAAAAAGAGCAAGAAAAGCAATGAGGAATGGCATATCCCCATTATTATTGA
TAATACCAACCTCCACGCCTGGGAAATGAAGCCCTATGCAGTCATGGCACTTAAAAATACTATGAAGTT
ATATTCCGAGAACCTGACACTCGCTGGAAATTAACGTTCAAGAGTTAGCAAGAAGAAACATTATGTTG
TCTCAAGAGAAAAATCCACCGAATGAAAGAACGGTATGAACACGATGTTACTTTTCACAGTGTGCTTCA
TGCAGAAAAGCCAAGCAGAATGAACAGAAACCAGGACAGGAATAATGCATTGCCTTCCAACAATGCCAGA
TACTGGAATTCCTACACAGAGTTTCCAAACCGGAGGGCCACGGTGGATTACAAATGAGAGCTCTATC
ACAGAAGGGCGGTTGTACCATGGATAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>RC224407 representing NM_052818

Red=Cloning site Green=Tags(s)

MEDSFLQSFGRSLQPPQQQQRQPRPPRGTTPRRHSFRKHLYLRLGLPGSGKTTLARQLQHDFPRAL
IFSTDDFFREDGAYEFNPDLFEEAHEWNQKRARKAMRNGISPIIIDNTNLHAWEMKPYAVMALENNYEV
IFREPDRWKFNVQELARRNIHGVSREKIHMKERYEHDVTFHSLHAEKPSRMNRNQDRNNALPSNNAR
YWNSYTEFPNRRAHGGFTNESSYHRRGGCHHG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms:

https://cdn.origene.com/chromatograms/mk6493_h07.zip

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_052818

ORF Size:

729 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_052818.3](#)

RefSeq Size: 3051 bp

RefSeq ORF: 732 bp


Locus ID: 90634

UniProt ID: [Q5TBK1](#)


Cytogenetics: 13q13.1

MW: 28.8 kDa


Product images:



Circular map for RC224407



Western blot validation of overexpression lysate (Cat# [LY409475]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from un-transfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC224407 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified N4BP2L1 protein (Cat# [TP324407]). The protein was produced from HEK293T cells transfected with N4BP2L1 cDNA clone (Cat# RC224407) using MegaTran 2.0 (Cat# [TT210002]).