

Product datasheet for **RC224412L1V**

BLOC1S1 (NM_001487) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BLOC1S1 (NM_001487) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BLOC1S1
Synonyms:	BLOS1; BORCS1; GCN5L1; MICoA; RT14
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001487
ORF Size:	375 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224412).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001487.1
RefSeq Size:	590 bp
RefSeq ORF:	462 bp
Locus ID:	2647
UniProt ID:	P78537
Cytogenetics:	12q13.2
MW:	14.3 kDa



[View online »](#)

Gene Summary:

BLOC1S1 is a component of the ubiquitously expressed BLOC1 multisubunit protein complex. BLOC1 is required for normal biogenesis of specialized organelles of the endosomal-lysosomal system, such as melanosomes and platelet dense granules (Starcevic and Dell'Angelica, 2004 [PubMed 15102850]).[supplied by OMIM, Mar 2008]