

Product datasheet for **RC224430L2V**

TULP4 (NM_020245) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	TULP4
Synonyms:	TUSP
Mammalian Cell	None
Selection:	
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_020245
ORF Size:	4629 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224430).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_020245.3 , NP_064630.2
RefSeq Size:	11123 bp
RefSeq ORF:	4632 bp
Locus ID:	56995
UniProt ID:	Q9NRJ4
Cytogenetics:	6q25.3



Protein Families:	Druggable Genome, Transcription Factors
MW:	168.8 kDa
Gene Summary:	May be a substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.[UniProtKB/Swiss-Prot Function]