

Product datasheet for **RC224430L3V**

TULP4 (NM_020245) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: TULP4

Synonyms: TUSP

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: NM_020245

ORF Size: 4629 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224430).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_020245.3](#), [NP_064630.2](#)

RefSeq Size: 11123 bp

RefSeq ORF: 4632 bp

Locus ID: 56995

UniProt ID: [Q9NRJ4](#)

Cytogenetics: 6q25.3



Protein Families: Druggable Genome, Transcription Factors

MW: 168.8 kDa

Gene Summary: May be a substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.[UniProtKB/Swiss-Prot Function]