

## Product datasheet for RC224434L2

### SR140 (U2SURP) (NM\_001080415) Human Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	SR140 (U2SURP) (NM_001080415) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	SR140
Synonyms:	fSAPa; SR140
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224434).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.

ACCN:	NM_001080415
ORF Size:	3087 bp



[View online »](#)

