

## Product datasheet for **RC224460L3V**

### eEF1A2 binding protein (IGFN1) (NM\_178275) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	eEF1A2 binding protein (IGFN1) (NM_178275) Human Tagged ORF Clone Lentiviral Particle
Symbol:	eEF1A2 binding protein
Synonyms:	DKFZp434B1231; eEF1A2 binding protein; EEF1A2BP1; immunoglobulin-like and fibronectin type III domain containing 1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_178275
ORF Size:	2604 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224460).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_178275.3</a>
RefSeq Size:	3951 bp
RefSeq ORF:	2606 bp
Locus ID:	91156
Cytogenetics:	1q32.1
MW:	94.1 kDa


[View online »](#)