

## Product datasheet for **RC224474L3V**

### **ARF1 (NM\_001024226) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | ARF1 (NM_001024226) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | ARF1   |
| Synonyms:                 | PVNH8  |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_001024226   |
| ORF Size:                 | 543 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC224474).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001024226.1</a> , <a href="#">NP_001019397.1</a>  |
| RefSeq Size:              | 1986 bp  |
| RefSeq ORF:               | 546 bp   |
| Locus ID:                 | 375  |
| UniProt ID:               | <a href="#">P84077</a>   |
| Cytogenetics:             | 1q42.13  |
| Protein Pathways:         | Vibrio cholerae infection  |
| MW:                       | 20.7 kDa   |



[View online »](#)

**Gene Summary:**

ADP-ribosylation factor 1 (ARF1) is a member of the human ARF gene family. The family members encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking as activators of phospholipase D. The gene products, including 6 ARF proteins and 11 ARF-like proteins, constitute a family of the RAS superfamily. The ARF proteins are categorized as class I (ARF1, ARF2 and ARF3), class II (ARF4 and ARF5) and class III (ARF6), and members of each class share a common gene organization. The ARF1 protein is localized to the Golgi apparatus and has a central role in intra-Golgi transport. Multiple alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]