

## Product datasheet for **RC224497L1V**

### VCX (NM\_013452) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                                  |                                      |
|----------------------------------|--------------------------------------|
| <b>Product Type:</b>             | Lentiviral Particles                 |
| <b>Symbol:</b>                   | VCX                                  |
| <b>Synonyms:</b>                 | VCX-10r; VCX-B1; VCX1; VCX10R; VCXB1 |
| <b>Mammalian Cell Selection:</b> | None                                 |
| <b>Vector:</b>                   | pLenti-C-Myc-DDK (PS100064)          |
| <b>Tag:</b>                      | Myc-DDK                              |
| <b>ACCN:</b>                     | NM_013452                            |
| <b>ORF Size:</b>                 | 618 bp                               |

**ORF Nucleotide Sequence:** The ORF insert of this clone is exactly the same as(RC224497).

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

|                      |                                                           |
|----------------------|-----------------------------------------------------------|
| <b>RefSeq:</b>       | <a href="#">NM_013452.2</a> , <a href="#">NP_038480.2</a> |
| <b>RefSeq Size:</b>  | 976 bp                                                    |
| <b>RefSeq ORF:</b>   | 621 bp                                                    |
| <b>Locus ID:</b>     | 26609                                                     |
| <b>UniProt ID:</b>   | <a href="#">Q9H320</a>                                    |
| <b>Cytogenetics:</b> | Xp22.31                                                   |



**MW:** 22.3 kDa

**Gene Summary:**

This gene belongs to the VCX/Y gene family, which has multiple members on both X and Y chromosomes, and all are expressed exclusively in male germ cells. The X-linked members are clustered on chromosome Xp22 and Y-linked members are two identical copies of the gene within a palindromic region on Yq11. The family members share a high degree of sequence identity, with the exception that a 30-bp unit is tandemly repeated in X-linked members but occurs only once in Y-linked members. The VCX gene cluster is polymorphic in terms of copy number; different individuals may have a different number of VCX genes. VCX/Y genes encode small and highly charged proteins of unknown function. The presence of a putative bipartite nuclear localization signal suggests that VCX/Y members are nuclear proteins. This gene contains 10 repeats of the 30-bp unit. [provided by RefSeq, Jul 2008]