

## Product datasheet for **RC224499L3V**

### **RABEP1 (NM\_001083585) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RABEP1 (NM_001083585) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RABEP1
Synonyms:	RAB5EP; RABPT5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001083585
ORF Size:	2487 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224499).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001083585.1</a>
RefSeq Size:	5322 bp
RefSeq ORF:	2490 bp
Locus ID:	9135
UniProt ID:	<a href="#">Q15276</a>
Cytogenetics:	17p13.2
Protein Pathways:	Endocytosis
MW:	95.3 kDa



[View online »](#)

**Gene Summary:**

Rab effector protein acting as linker between gamma-adaptin, RAB4A and RAB5A. Involved in endocytic membrane fusion and membrane trafficking of recycling endosomes. Involved in KCNH1 channels trafficking to and from the cell membrane (PubMed:22841712). Stimulates RABGEF1 mediated nucleotide exchange on RAB5A. Mediates the traffic of PKD1:PKD2 complex from the endoplasmic reticulum through the Golgi to the cilium (By similarity). [UniProtKB/Swiss-Prot Function]