

Product datasheet for RC224500

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com

OriGene Technologies, Inc.

techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Caspase 10 (CASP10) (NM_032974) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Caspase 10 (CASP10) (NM_032974) Human Tagged ORF Clone

Tag: Myc-DDK
Symbol: Caspase 10

Synonyms: ALPS2; FLICE-2; FLICE2; MCH4

Mammalian Cell Neomycin

Selection:


Vector:pCMV6-Entry (PS100001)E. coli Selection:Kanamycin (25 ug/mL)

Restriction Sites: Sgfl-Mlul

Cloning Scheme:

Cloning sites used for ORF Shuttling:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_032974

ORF Size: 1563 bp



Caspase 10 (CASP10) (NM_032974) Human Tagged ORF Clone - RC224500

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 032974.5</u>

RefSeq Size:2054 bpRefSeq ORF:1566 bp

 UniProt ID:
 Q92851

 Cytogenetics:
 2q33.1

Locus ID:

Protein Families: Druggable Genome, Protease

843

Protein Pathways: Apoptosis, RIG-I-like receptor signaling pathway

MW: 58.8 kDa

Gene Summary: This gene encodes a protein which is a member of the cysteine-aspartic acid protease

(caspase) family. Sequential activation of caspases plays a central role in the execution-phase


of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein cleaves and activates caspases 3 and 7, and the protein itself is processed by caspase 8. Mutations in this gene are associated with type IIA autoimmune lymphoproliferative syndrome, non-Hodgkin lymphoma and gastric cancer.

Alternatively spliced transcript variants encoding different isoforms have been described for

this gene. [provided by RefSeq, Apr 2011]



Product images:



Circular map for RC224500