

Product datasheet for **RC224505L4V**

CD200RIL (NM_001008784) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	CD200RIL
Synonyms:	CD200R2; CD200RLa
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001008784
ORF Size:	813 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224505).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001008784.1
RefSeq Size:	1296 bp
RefSeq ORF:	816 bp
Locus ID:	344807
UniProt ID:	Q6Q8B3
Cytogenetics:	3q13.2



Protein Families: Druggable Genome, Transmembrane

MW: 29.9 kDa

Gene Summary: May be a receptor for the CD200/OX2 cell surface glycoprotein.[UniProtKB/Swiss-Prot Function]