

Product datasheet for **RC224513L4V**

MOBKL2C (MOB3C) (NM_201403) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MOBKL2C (MOB3C) (NM_201403) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MOBKL2C
Synonyms:	MOB1E; MOBKL2C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_201403
ORF Size:	648 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC224513).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_201403.1
RefSeq Size:	2822 bp
RefSeq ORF:	651 bp
Locus ID:	148932
UniProt ID:	Q70IA8
Cytogenetics:	1p33
MW:	25.6 kDa


[View online »](#)

Gene Summary:

The protein encoded by this gene is similar to the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. Alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]