

## Product datasheet for **RC224514A1V**

### Human Saitohin (STH) (NM\_001007532) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human Saitohin (STH) (NM\_001007532) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Saitohin  
**Synonyms:** MAPTIT  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC224514 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAGTGAGGGTGGAGGCCAAGTCTCATGCATTTTTGCAGCCCCACAAGACTGTGCAGGTGGCCGGCC  
TCATTGAATGCGGGTTAATTTAACTCAGCCTCTGTGTGAGTGGATGATTCAGGTTGCCAGAGACAGAAC  
CCTCAGCTTAGCATGGGAAGTAGCTTCCCTGTTGACCCGTAGTTCATCTGAGGTTGGCTTGAAGGTGTG  
GGCACCATTTGGCCAGTTCTTACAGCTCTGAAGAGAGCAGCAGGAATGGGCTGAGCAGGGAAGACAAC  
TTTCCATTGAAGGCCCTTTAGGGCCAGAAGTCCCTCCCACCTGCAGCTGCCCTGCCTCTGCCCAT  
GAGGGTGAGAGTCAGGCGACCTCATGCCAAGTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC224514 protein sequence  
Red=Cloning site Green=Tags(s)

MSEGGQVSCIFAAPTRLCRWPALIECGVNLTPQLCEWMIQVARDRTLALWEVASLLTSSSEVGLGV  
GTIWPSSYSSESSRNQAEQGRQLSIEGPFQGNCPSPHAAALPLPMRGESQATSCQV

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001007532



[View online »](#)

|                      |   |
|----------------------|---|
| <b>ORF Size:</b>     | 384 bp  |
| <b>Buffer:</b>       | PBS with 0.001% Pluronic F68  |
| <b>Stability:</b>    | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| <b>RefSeq:</b>       | <u><a href="#">NM_001007532.1</a></u>   |
| <b>RefSeq Size:</b>  | 445 bp  |
| <b>RefSeq ORF:</b>   | 387 bp  |
| <b>Locus ID:</b>     | 246744  |
| <b>UniProt ID:</b>   | <u><a href="#">Q8IWL8</a></u>   |
| <b>Cytogenetics:</b> | 17q21.31  |
| <b>MW:</b>           | 13.7 kDa  |