

## Product datasheet for **RC224525A1V**

### Human Gemin 7 (GEMIN7) (NM\_001007270) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human Gemin 7 (GEMIN7) (NM\_001007270) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Gemin 7  
**Synonyms:** SIP3  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC224525 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGCAAACTCCAGTGAACATTCCCGTGCCTGTGCTCCGGCTGCCCCGGGGCCCTGATGGCTTCAGCCGTG  
 GCTTTGCCCTGATGGACGCAGAGCCCCCTTGAGGCCAGAGGTTCTGAAATCCAGGAGTGTCCCATAGC  
 TCAAGAATCCCTGGAATCCCAGGAGCAGCGGGCAGAGCCGCCCTTCGGGAGCGTTACCTCCGCAGCCTG  
 CTGGCCATGGTGGGTATCAGGTGAGCTTACGTTGCACGAGGGTGTGCGTGTGGCCGCCCACTTTGGAG  
 CCACCGACCTGGATGTGGCCAACTTCTACGTGTCACAGCTGCAGACTCCCATAGGTGTGAAGCAGAGGC  
 GCTGCTCCGATGTAGTGACATTATTTTCATATACCTTCAAGCCA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC224525 protein sequence  
 Red=Cloning site Green=Tags(s)

MQTPVNI PVPVLR LPRGPDGFSR GFAPDGRRAPLRPEVPEIQECP IAQESLESQEQRARAALRERYLRSL  
 LAMVGHQVSFTLHEGVVAAHFGATDL DVANFYVSQLQTPIGVQAEALLRCSDIISYTFKP

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001007270



[View online »](#)

<b>ORF Size:</b>	393 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_001007270.1</a></u>
<b>RefSeq Size:</b>	1561 bp
<b>RefSeq ORF:</b>	396 bp
<b>Locus ID:</b>	79760
<b>UniProt ID:</b>	<u><a href="#">Q9H840</a></u>
<b>Cytogenetics:</b>	19q13.32
<b>MW:</b>	14.5 kDa