

Product datasheet for **RC224536L2V**

VPS18 (NM_020857) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	VPS18
Synonyms:	PEP3
Mammalian Cell	None
Selection:	
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_020857
ORF Size:	2919 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC224536).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_020857.2
RefSeq Size:	3922 bp
RefSeq ORF:	2922 bp
Locus ID:	57617
UniProt ID:	Q9P253
Cytogenetics:	15q15.1



Domains: Pep3_Vps18

MW: 110 kDa

Gene Summary: Vesicle mediated protein sorting plays an important role in segregation of intracellular molecules into distinct organelles. Genetic studies in yeast have identified more than 40 vacuolar protein sorting (VPS) genes involved in vesicle transport to vacuoles. This gene encodes the human homolog of yeast class C Vps18 protein. The mammalian class C Vps proteins are predominantly associated with late endosomes/lysosomes, and like their yeast counterparts, may mediate vesicle trafficking steps in the endosome/lysosome pathway. [provided by RefSeq, Jul 2008]